PROTECTIVE MASK

GB2626-2006 KN95

EU 149:2001 +A1:2009 • FFP2 NR D

FFP2

PROTECTIVE MASK

GB2626-2006 (Civil Grade)

KN95

APPROVED

Antiviral fiber

Skin-friendly

Particulate matter
filtration

Individual package

FDA

CE

EU 149:2001 +A1:2009 • FFP2 NR D

Product Parameters

Standard size/Earloop

**Individual Package
Clean & Hygienic**

Product Name: PRD Protective Masks

Main Material: Nonwoven Fabric, Meltblown Fabric,
Needle-punched Cotton

Execution Standard: GB2626-2006 / EU 149:2001 +A1:2009

Expiration Date: 3 Years

See the packaging indication for production batch number
and expiration date.

Manufacturer: DONGGUAN LEIHUO MEDICAL DEVICE CO.,LTD

Address: 2F, Building B & 3F, Building D, No.30 Zhonghe Industrial
Park, Liusha, Puning city, Jieyan city, Guangdong, Province

Place of Production: Dongguan Guangdong

PROTECTION LEVEL

KN95

FFP2

N95

EU 149:2001 +A1:2009 · FFP2 NR D

Product Parameters


SPECIFICATIONS

MASK

(Independent Packages)

Size:10.7x16cm


40 pcs/Box

Size:16x13x23cm

MEANS:L*W*H
 69*37*27CM
 G.W: 5.3KG
 N.W:3.2KG
 1 CASE: 10 boxes
 1 BOX: 40pcs
 1 CASE: 400pcs

New & Good Materials, No Secondary Pollution

Multi-effect Protection & Filtration

Product Parameters

Product Name: PRD Protective Masks

Model: LH-KN95(KHA)

Main Material: Nonwoven Fabric, Meltblown Fabric, Needle-punched Cotton
Protection Level: N95

Execution Standard: GB2626-2006 / EU 149: 2001 +A1: 2009

Expiration Date: 3 Years

See the packaging indication for production batch number and expiration date.

Manufacturer: DONGGUAN LEIHUO MEDICAL DEVICE CO.,LTD

Address: Dongshan Yongshenglu No.47 Room 101, Qishi Dongguan,
Guangdong, 523500, China

Place of Production: Dongguan Guangdong

Precautions

1. If the package is damaged, it is forbidden to use the mask. It is limited to respiratory protection against non-oily particles only.

2. It is recommended to store in a dry, ventilated, non-corrosive gas environment below a humidity of 80%, Keep away from ignition sources and flammable materials.

3. This product is not applicable to special industries such as anoxic environment, underwater operation, fire control and industrial dust protection or to respiratory protection of infants.

4. Do not wear a mask if unventilated, hard in breathing or during sleep.

5. People with abnormal heart and lung functions should use masks with caution.

6. If a mask is damaged or you feel the respiratory resistance is obviously increased, replace it in time.

7. In case of any discomfort or adverse reactions during use, discontinue use immediately.

8. Do not reuse.

Nonapplicable Scenes

- Nonapplicable Scenes
- Do not use in an environment above 50°C.
- Do not use in an environment with an oxygen content of less than 20%.
- Do not use in a toxic gas environment.
- Not recommended for children, pregnant women or the elderly.

MADE IN C.C.T

How to Use

				
<p>1. Open the mask to make the nose clip at the top, and pull the ear straps with both your hands.</p> <p>2. Hold the mask against your chin to completely cover your nose and mouth.</p> <p>3. Pull the ear straps behind your ears and adjust them to make you feel comfortable.</p>			<p>4. Use both your hands to adjust the shape of the nose clip. Place your fingers in the middle of the nose clip and press it inwards while moving your fingertips along both sides of the nose clip until it is pressed to fit the bridge of your nose. (Making the nose clip get a seal with only one hand may affect the tightness of the mask).</p> <p>5. Cover the mask with your hand and exhale vigorously. If you feel the air escaping from the nose clip, it is required to tighten the nose clip; if the air escapes from the edge of the mask, readjust the headband to ensure tightness.</p>	

Inspection Method

Before use, check whether the package of the mask is damaged and the components are intact. In case of any stain or damage, discontinue use immediately. Before entering the work area, it is required to check the air tightness of the mask. If you feel a leak, please wear it again according to the instructions for how to wear a mask. You can enter the work area only after feeling no leak. (If you feel the respiratory resistance is obviously increased or the mask falls short of the standard for use, it is recommended to replace the mask immediately.)

KN95	GB2626-2006 (Civil Grade)	This dust mask is used to prevent non-oily suspended particles. Filtration Efficiency $\geq 95\%$.
FFP2 N95		

Do not use if package is damaged!
Please refer to the manual before
use and use a mask correctly!

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Brand Name : Leihuo
Type : Protective masks
Model N : LH-KN95

Manufacturer : Dongguan Leihuo Medical Device Co., Ltd.
Address : Dongshan Yongshenglu N0.47 Room 101, Qishi, Dongguan,
Guangdong China

Standard Used : EN 14683:2005

Conclusion :

After inspection of the technical documentation issued by the customer, and in his request, we express our opinion that the product meets the technical requirement of the following directives and standards:

93/42/EEC Medical devices (MDD)

This opinion is only valid for the directive, the equipment and configuration described, in conjunction with the test data detailed above and with compliance with all applicable legal requirement for the product .

The following manufacturer documents was inspected:

Presence of Declaration of conformity template	✓ OK
Presence of test report using standards as indicated in the declaration of conformity Test report reference : HTT202003147LR	✓ OK
Presence of  symbol in the product label.	✓ OK
Presence of instruction manual	✓ OK
Use of valid Harmonized standard in the declaration of conformity	✓ OK
Presence of product description in the technical construction file	✓ OK

Copyright of this Certificate is owned by CELAB® Italy and may not be reproduced other than in full and with the prior approval of the General Manager. Use of this certificate is subjected to Celab regulation available on Celab web site.

Check the authenticity of this certificate and related information before use in the web site www.celab.com introducing the UCN number in the 'Check document authenticity' area. You will see copy of this certificate and regulation on certificate use. This document is released only for scope allowed by laws- Do not use this document without full understanding of regulation.

Massimiliano Bertoldi
General Manager – CELAB
www.celab.com

www.celab.com

Fiscal Year 2020

CERTIFICATION OF FDA REGISTRATION

This certifies that:

DONGGUAN LEIHUO MEDICAL DEVICE CO., LTD.
Dongshan Yongshenglu No.47 Room 101, Qishi
Dongguan, Guangdong, 523500, CHINA

has completed the FDA Establishment Registration and Device Listing with the US Food & Drug Administration, through

Shenzhen CCT Testing Technology Co., Ltd.

Owner/Operator Number: 10063202

Device Listing#:


Listing No	Code	Device Name	Proprietary Name
D375637	KHA	MASK, SCAVENGING	KN95 Mask, Disposable Mask

CCT will confirm that such registration remains effective upon request and presentation of this certificate until the end of the calendar year stated above, unless said registration is terminated after issuance of this certificate. CCT makes no other representations or warranties, nor does this certificate make any representations or warranties to any person or entity other than the named certificate holder, for whose sole benefit it is issued. This certificate does not denote endorsement or approval of the certificate-holder's device or establishment by the U.S. Food and Drug Administration. CCT assumes no liability to any person or entity in connection with the foregoing.

Pursuant to 21 CFR 807.39, "Registration of a device establishment or assignment of a registration number does not in any way denote approval of the establishment or its products. Any representation that creates an impression of official approval because of registration or possession of a registration number is misleading and constitutes misbranding." The U.S. Food and Drug Administration does not issue a certificate of registration, nor does the U.S. Food and Drug Administration recognize a certificate of registration, CCT is not affiliated with the U.S. Food and Drug Administration.

Shenzhen CCT Testing Technology Co., Ltd.
W: www.fda-test.com E: fda@fda-test.com
T: 400-8788-298 T: 86-755-36916737

Web: <http://www.fda.gov> Tel: 1-888-INFO-FDA (1-888-463-6332) e-mail: webmail@oc.fda.gov

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